

## CHIRASCAN AND CHIRASCAN V100 VERSATILE RESEARCH SYSTEMS




### HIGH PERFORMANCE, READY TO RUN

- Determine structural and thermodynamic properties
  - Gain insight and detect changes in secondary and tertiary structure
  - Determine response to thermal or chemical changes
  - Study folding and unfolding mechanisms
- Achieve highest sensitivity and accuracy
- Generate highest quality data
- Optimize sample concentration and absorbance
- Expand capabilities with dedicated Chirascan accessories

# DETERMINE STRUCTURAL AND THERMODYNAMIC PROPERTIES


## Gain insight and detect changes in secondary and tertiary structure

Secondary structure: far-UV spectrum of a globular protein



Simultaneous acquisition of CD and absorbance spectra, 0.5 mm pathlength, Chirascan V100. Courtesy of leading research university, Germany

Tertiary structure: near-UV spectra of two monoclonal antibodies



Differences between near-UV spectra due to slight changes in orientation of aromatic moieties, Chirascan V100, 10 mm pathlength



## Determine thermodynamic properties - continuous thermal ramping

- Monitor at each wavelength
- Typical run: 70 spectra in 70 min, 1°C/min
- Record temperature directly - thermocouple in sample
- Derive melting points and enthalpies for multiple thermal transitions
- Associate change in structure with each thermal transition

pH	Melting temperature (°C)	van't Hoff enthalpy (kJ/mol)
pH 2	55.4	354
pH 3	69.4	385
pH 4	75.8	380
pH 5	76.9	400
pH 6	74.2	423
pH 7	72.7	367

Six datasets analyzed using Chirascan global thermodynamic analysis

Effect of pH on thermal denaturation



Two of six denaturation datasets acquired at pH 2-7, lysozyme, Chirascan 6-cell turret, Chirascan V100, raw data, no baseline adjustment, no smoothing, 0.5 mm pathlength


# ACHIEVE HIGHEST SENSITIVITY AND ACCURACY

Since their introduction in 2005, Chirascan™ systems have continued to feature in thousands of peer-reviewed publications covering a wide range of research areas. Chirascan V100 now offers the increased sensitivity and accuracy preferred for CD analysis of biomolecules.

- Avalanche photodiode detector enhances sensitivity
- Increased signal:noise compared to conventional photomultiplier
- Accurate normalization from simultaneous measurement of absorbance and CD

Increased sensitivity when sample is limited


Chirascan V100 (avalanche photodiode detector)  
Chirascan (photomultiplier)



Tertiary structure of lysozyme – raw data, no smoothing, 10 min. baseline / 10 min. sampling, n=3 scans, 0.5 nm step, 10 mm pathlength, spectra offset for clarity


Increased sensitivity: faster measurements for thermal studies or photolabile samples

Chirascan V100 (avalanche photodiode detector)  
Chirascan (photomultiplier)



Tertiary structure of lysozyme – raw data, no smoothing, baseline corrected, n=3 scans, 1 nm step, 10 mm pathlength, spectra offset for clarity

- Accurate CD values across entire wavelength range
- Overcome challenges of chemical calibration
- Optics-based, multiwavelength calibration



Conventional chemical calibration methods require considerable skill in preparation. Standards, such as camphor-10-sulfonic acid (CSA), are unstable, photolabile and hygroscopic. In addition, single wavelength calibration (290.5 nm) assumes the same linear response at all wavelengths.

The optics-based, multiwavelength calibration method used in Chirascan V100 overcomes these challenges. The correct calibration is applied to every wavelength to yield accurate CD values.

## READY TO RUN – GENERATE HIGHEST QUALITY DATA

Chirascan systems are supplied with features and accessories required for acquisition of high quality CD data – from built-in temperature control during analysis to cuvettes for the most common analytical conditions.\* A basic training program follows installation to familiarize users new to Chirascan.

### PHOTOELASTIC MODULATOR

- Converts horizontally polarized light to circularly polarized light. Alternates between left- and right-handed circular polarized light

### MONOCHROMATOR

- Produces horizontally, linearly polarized monochromatic light
- Two polarizing prisms maximize light throughput

### AIR-COOLED XENON LAMP

- Software-controlled
- Up-time recorded

### ACTIVE NITROGEN MANAGEMENT SYSTEM

- Regulates purge gas consumption
- Software-controlled

### AVALANCHE PHOTODIODE DETECTOR (CHIRASCAN V100)

- Highest sensitivity (high signal: noise)

### PHOTOMULTIPLIER DETECTOR (CHIRASCAN)



### MOLECULAR SIEVE, ACTIVATED CHARCOAL FILTER

- Removes common gas impurities

\* Your local Applied Photophysics representative can supply specific details of components supplied for your region.



## TEMPERATURE-CONTROLLED SAMPLE CHAMBER

- Consistent analytical conditions
- Continuous temperature ramps
- Temperature measured directly in sample

## OPTICS-BASED, MULTIWAVELENGTH CALIBRATION

- For CD accuracy (Chirascan V100)

## CUVETTES AND HOLDERS


- Selected for far- and near-UV CD analysis of biomolecules



## WATER CIRCULATOR

- Dissipates heat from sample chamber Peltier

# CONTROL AND ANALYSIS SOFTWARE



## CHIRASCAN CONTROL

- Easily define run parameters and store routine protocols
- Saves time with scheduled start-up/shutdown of lamp and  $N_2$  supply
- Fail-safe lamp switch-off if  $N_2$  flow drops
- Ensures  $O_2$ -free conditions with  $N_2$  purge

## GLOBAL THERMODYNAMIC ANALYSIS

- Derive melting points and enthalpies from multiwavelength, thermal denaturation experiments

## OPTIMIZE SAMPLE CONCENTRATION AND ABSORBANCE: CUVETTES AND HOLDERS

Selecting a suitable cuvette with optimal pathlength is critical to acquisition of highest quality data. Cuvettes for Chirascan systems are manufactured from far-UV quartz to enable analysis of secondary structure. The range of cuvettes and compatible holders provides full flexibility when optimizing sample concentration and absorbance.

For secondary structure (far-UV) analysis	
0.5 mm 175 $\mu$ L one-piece stoppered cuvette Adaptor for 0.5 mm and 1 mm one-piece cuvettes	Supplied with Chirascan and Chirascan V100 Not suitable for fluorescence
1.0 mm 350 $\mu$ L one-piece stoppered cuvette	Not suitable for fluorescence; requires adaptor
2.0 mm 700 $\mu$ L one-piece stoppered cuvette Spacer for 2 mm pathlength cuvette	Not suitable for fluorescence
2.0 mm 400 $\mu$ L one-piece stoppered cuvette for manual systems	Simultaneous measurement of secondary structure by CD and tertiary structure by fluorescence. No adaptor needed. Not compatible with Chirascan 6-cell turret.
4.0 mm 1400 $\mu$ L one-piece stoppered cuvette	Simultaneous measurement of secondary structure by CD and tertiary structure by fluorescence. No adaptor needed.
0.1 mm 160 $\mu$ L cylindrical cell, Holder for cylindrical cells	Useful for highly absorbing/chiral buffers. Scan further into far-UV at fixed pathlength. Not suitable for fluorescence.
For tertiary structure (near-UV) analysis	
10 mm 3500 $\mu$ L one-piece stoppered cuvette	Supplied with Chirascan and Chirascan V100 Suitable for fluorescence, certified free from strain birefringence.
10 mm 600 $\mu$ L one-piece stoppered cuvette Holder for 10 mm 600 $\mu$ L cuvettes	Suitable for fluorescence; reduced volume
5.0 mm 1750 $\mu$ L one-piece stoppered cuvette Spacer for 5 mm pathlength cuvettes	Not suitable for fluorescence.

For secondary structure analysis with rapid cell cleaning. Scan further into the far-UV Not suitable for fluorescence or for use with chiral buffers	
Adaptor for demountable/slide cells	
0.01 mm 3 $\mu$ L demountable/slide cell	
0.1 mm 30 $\mu$ L demountable/slide cell	
0.2 mm 60 $\mu$ L demountable/slide cell	
0.5 mm 150 $\mu$ L demountable/slide cell	

# EXPAND CAPABILITIES WITH DEDICATED CHIRASCAN ACCESSORIES

## CCD FLUOROMETER

Monitor changes in fluorescence



## 6-CELL TURRET

Increase capacity and productivity



## TITRATOR AND PH PROBE

Monitor concentration- and pH-dependent changes in CD, fluorescence or absorbance



## LD COUETTE CELL

Use linear dichroism to gain insight into conformation and orientation of molecular structures



## STOPPED-FLOW

Characterize fast reactions, complement CD spectra with kinetic information



## INTEGRATING SPHERE AND SOLID SAMPLE HOLDER

See changes in chirality of solid state samples



## ALSO AVAILABLE

A fluorescence anisotropy detector to characterize ligand binding events.  
An optical rotatory dispersion accessory to characterize chiral molecules.  
A magnetic CD accessory for study of molecules such as metalloproteins.  
A near IR extension kit to expand scanning range of a Chirascan V100 to 1700 nm.

Please contact your Applied Photophysics representative to discuss requirements, and system compatibility.